

NOTHING CONNECTS QUITE LIKE THE ORIGINAL!

SECUR-LOK* WITH SAFETY RING

- Excellent fit into the MIC-KEY* feeding port, opening or unlocking the anti-reflux valve
- Less risk of leakage due to the designed and tested compatibility
- Safety ring helps prevent inadvertent connection to intravenous lines



A COMPREHENSIVE RANGE OF MIC-KEY* EXTENSION SETS

The SECUR-LOK* connector with safety ring at the patient end fits perfectly into the MIC-KEY* feeding port without damaging the anti-reflux valve.

FOR CONTINUOUS FEEDING:

Products 0141-12 and 0141-24

- Featuring male ENFit® feeding and medication connectors with screw caps
- Available in 2 lengths: 30 and 60 cm

Products 0144-12 and 0144-24

- Equipped with a single ended male ENFit® feeding connector
- Available in 2 lengths: 30 and 60 cm

FOR BOLUS FEEDING: products 0143-12 and 0143-24

- Wide bore tubing allows hydration and bolus feeding with a Female ENFit® syringe
- Available in 2 lengths: 30 and 60 cm

FOR ADMINISTRATION OF MEDICATION: product 0142-02

- The 5 cm long extension set with double ENFit® connector and screw cap is especially useful to administer medication

The SECUR-LOK* connector with safety ring at the patient ends fits perfectly into the MIC-KEY* feeding port without damaging the anti-reflux valve. All enteral feed, medication and fluids must be administered via a MIC-KEY* Extension Set connected to the MIC-KEY* feeding tube.



MIC-KEY* G -TUBE: MAKE SURE YOU OBTAIN THE ORIGINAL!

ASK FOR THE MIC-KEY* EXTENSION SET
WITH **SECUR-LOK*** AND **SAFETY RING.**



Wash the extension set after every feeding with warm soapy water and rinse it thoroughly. Store in a dry and clean place.

Prompt flushing and rinsing prevents the tube feeds from drying and building up. For optimum functionality it is recommended that you replace your MIC-KEY* Extension Sets every 2 weeks.



DESCRIPTION	LENGTH	PACKAGING	PRODUCT CODE
MIC-KEY* Extension Set ENFit® with Right-Angle Connector	30 cm	5 Eaches	0141-12
MIC-KEY* Extension Set ENFit® with Right-Angle Connector	60 cm	5 Eaches	0141-24
MIC-KEY* Medication Set ENFit® with Right-Angle Connector	5 cm	5 Eaches	0142-02
MIC-KEY* Bolus Extension Set ENFit® with Straight Connector	30 cm	5 Eaches	0143-12
MIC-KEY* Bolus Extension Set ENFit® with Straight Connector	60 cm	5 Eaches	0143-24
MIC-KEY* Extension Set single ended ENFit® with Right-Angle Connector	30 cm	5 Eaches	0144-12
MIC-KEY* Extension Set single ended ENFit® with Right-Angle Connector	60 cm	5 Eaches	0144-24

There are inherent risks in all medical devices, please refer to the product labeling for Indications, Cautions, Warnings and Contraindications.

AVANOS

For more information, please send an email to customerservice.ANZ@avanos.com or visit www.avanos.com.au.

Alternatively, more information and resources can be found at tubefed.com.au

* Registered Trademark or Trademark of Avanos Medical, Inc., or its affiliates. © 2018 AVNS. All rights reserved.
ENFit® is a Registered Trademark of Global Enteral Device Supplier Association, Inc. Used with their permission. COPY-05969